ClearBand 100x Phosphatase Inhibitor Cocktail

1 ml

Cat No: PHIC1

Shipping : Ship with blue ice. **Storage** : Store at -20°C.

General Information

ClearBand 100x Phosphatase Inhibitor Cocktail is formulated for efficient inhibition of endogenous phosphatases during cell lysis and protein extraction. *ClearBand* 100x Phosphatase Inhibitor Cocktail is a blend of 5 phosphatase inhibitors that can broadly inhibit acid phosphatases, alkaline phosphatases, and protein tyrosine phosphatases (PTPs).

Suitable for use in following applications;

WB, Co-IP, pull-down, IF, IHC, Flow Cytometry, kinase assay and etc.

Protocol

- 1. Thaw ClearBand 100x Phosphatase Inhibitor Cocktail on ice,
- 2. Add at 1:100 (v/v) dilution to solution samples (such as cell lysates or tissue extracts) before assaying.